



Important Safety Instructions:

1. Read all instructions before using the electric toothbrush.
2. Keep this manual for future reference.
3. Do not use the electric toothbrush if the cord or plug is damaged.
4. Do not immerse the electric toothbrush in water or any other liquid.
5. Do not use the electric toothbrush if it is not working properly or has been dropped or damaged.
6. Do not use the electric toothbrush while bathing or showering.
7. Keep the electric toothbrush away from children and pets.
8. Do not use attachments or accessories that are not recommended by the manufacturer.
9. If the electric toothbrush is used by more than one person, replace the brush head after each use.

Charging the Electric Toothbrush:

1. Before first use, charge the electric toothbrush for at least 24 hours.
2. Connect the charging base to a power outlet.
3. Place the electric toothbrush on the charging base and ensure it is properly aligned.
4. The charging indicator light will illuminate to indicate that the electric toothbrush is charging.
5. Once fully charged, the charging indicator light will turn off.

Using the Electric Toothbrush:

1. Apply toothpaste to the brush head.
2. Wet the brush head and turn on the electric toothbrush by pressing the power button.
3. Gently move the brush head in a circular motion along the teeth and gums.
4. Do not press too hard as excessive pressure may damage the teeth and gums.
5. Brush for at least two minutes, spending equal time on each quadrant of the mouth.
6. Rinse the brush head and handle thoroughly after each use.
7. Store the electric toothbrush in a clean, dry place.

Replacing the Brush Head:

1. Replace the brush head every three months or sooner if the bristles become frayed or worn.

2. To remove the brush head, firmly grip the base and pull the brush head away from the handle.
3. Align the new brush head with the handle and push it firmly until it clicks into place.

Troubleshooting:

1. If the electric toothbrush does not turn on, ensure it is properly charged.
2. If the electric toothbrush stops working during use, check if the battery is depleted and needs to be recharged.
3. If the electric toothbrush exhibits abnormal behavior or does not function properly, contact customer support for assistance.

Note: This electric toothbrush is intended for personal use only. Do not share the brush head with others.